

## FERROCHELATASE (FECH)

### Sample preparation.

Suspend ~50- $\mu$ L mammalian cell pellet in 400 $\mu$ L TGC buffer (Tris buffered glycerol with dithiothreitol (DTT), made by dissolving 2mL or 2.52 g glycerol and 1.5mg DTT in 8.0 mL 20mM Tris pH 8.0).

Sonicate (homogenize) while in an ice bath at the lowest practicable power setting for 3 cycles x 5 sec at 50% duty (2.5 sec on, 2.5 sec off).

Determine the protein content and dilute to 1 $\mu$ g protein/ $\mu$ L with more TGD.

### Ferrochelatase reaction.

Prepare three 50- $\mu$ L aliquots of the cell preparation, two live and one inactivated in boiling water for 10 minutes (as control for non-enzymatic product formation).

Prepare the incubation buffer containing 160mM Tris pH 8.0, 40mM Bicine pH 8.0, 10mg/ml Tween20 and 0.38mg/mL palmitic acid),

Mix 150 $\mu$ L incubation buffer and 25 $\mu$ L zinc substrate (1mM aqueous Zn acetate) and pre-incubate for 5 minutes at 37°C.

Mix each 50- $\mu$ L aliquot of cell preparation with the 150 $\mu$ L incubation buffer plus zinc substrate.

Then add 25  $\mu$ L of mesoporphyrin IX substrate (250 $\mu$ M in 160mM Tris pH 8.0, 40mM Bicine pH8.0, 2mg/ml Tween20).

Incubate the mixture for 30 min at 37°C.

Add 750 $\mu$ L stop reagent (270 $\mu$ M ethylenediaminetetraacetic acid in a mixture containing dimethylsulfoxide-methanol, 30/70 by volume respectively).

Cool on ice for 15-20 min.

Centrifuged at 1500xg for 10 min at room temperature.

### Quantitation of product Zn mesoporphyrin IX.

Inject 10 $\mu$ L of supernatant solution of porphyrins in 1.5M HCl into a Waters Acquity UPLC (ultra performance liquid chromatography ) system, which includes a binary solvent manager, sample manager, fluorescence detector (FLR), column heater and an Acquity UPLC BEH C18, 1.7  $\mu$ M, 2.1 x 100 mm column. Set the FLR for zinc mesoporphyrin IX (ZnMeso) at 406 nm excitation and 578 nm emission. Quantify the ferrochelatase product relative to a standard ZnMeso solution, also in the stop reagent. Keep the sample chamber dark and at

ambient temperature. Solvent A is 0.2% aqueous formic acid while Solvent B was 0.2% formic acid in methanol. Set the flow rate at 0.40 mL per minute at 60°C for the total run time of 7 min. Use the following successive gradient settings for run time in minutes versus A: 0.0, 80%; 2.5, 1%; 4.5, 1%; 5.0, 80%. Keep all solvent gradients are linear.