

Protocol for reverse FECH assay. Hector Bergonia 03-26-20

Sonicate the washed pellet of cultured bacterial cells in two volumes of 10mM KPi pH 5.5. Determine protein content and adjust to about 20mg/mL with the same buffer to prepare the homogenate. Deoxygenate all samples and reagents.

Mix about 50 μ L homogenate (1.0 mg total protein) with 50 μ L of assay reagent containing 100 μ M hemin-imidazole, 4 mM ascorbic acid in 10 mM potassium phosphate buffer pH 5.5 under argon atmosphere.

Cap the sample tubes securely, incubate for 1h at 45°C in a water bath and then add 400 μ L 50% v/v acetone in ethanol.

Centrifuge the samples at 13,500 g for 10 min.

Analyze the supernatant for PPIX fluorescence by UPLC (ultra-performance liquid chromatography).

For blank controls samples of bacterial homogenates were heated for 10 min in a boiling water bath and assayed with the live samples in pairs.